



# CURA F

## VERSATILE & REUSABLE

This reusable and user-friendly faceshield was specifically designed and developed for healthcare professionals. Compatible with all prescription glasses, CURA-F is perfectly suited to provide effective protection against viruses, droplets and splashes. Easily assembled, adjusted, and cleaned, it ensures long-lasting comfort, performance and reliability. CURA-F is designed, manufactured and assembled in the UK using locally sourced sustainable materials. Because protecting the planet means protecting lives in the long term.

+ **LIGHTWEIGHT DISINFECTABLE FRAME (POLYAMIDE)**

+ **TINTED, EASY PEEL-OFF OUTER PROTECTIVE FILM**

+ **WASHABLE AND REUSABLE**

+ **PLATINUM LITE, ANTI-FOG COATING**



### TECHNOLOGIES

PLATINUM LITE

### VERSIONS



PFSCURFP02



PFSCURFP03



PFSCURFP04



## CURA F

VERSION	REFERENCE	LENS MARKING	FRAME MARKING	COATING	STANDARD
CLEAR	PFSCURFP02	CE 38 Z87 CROWN	CE 38 Z87 D3	PLATINUM LITE	CE-UKCA
CLEAR	PFSCURFP03	CE 38 Z87 CROWN	CE 38 Z87 D3	PLATINUM LITE	CE-UKCA
CLEAR	PFSCURFP04	CE 38 Z87 CROWN	CE 38 Z87 D3	PLATINUM LITE	CE-UKCA

BOLLÉ SAFETY PRODUCT APPROVAL CERTIFICATES TO CE STANDARDS WILL BE SENT ON REQUEST.



PFSCURFP02

HEADBAND/TEMPLES/BUCKLES



PFSCURFP03

HEADBAND/TEMPLES




PFSCURFP04

HEADBAND/BUCKLES

## GENERAL SPECIFICATIONS

LENS	PET	STRAP	ISOPRENE - WHITE
FRAME	PA6 - HEADBAND/TEMPLES/BUCKLES	USE	MULTI USE
WEIGHT	68 g	SUPPLIED WITH	

KNOW MORE  
ABOUT  
MARKINGS:



OR SCAN:

## PACKAGING &amp; EAN13 CODE

SINGLE:	INNER BOX:	OUTER PACK:
Size: 38 x 27 x 2 mm	Size: xx cm Weight: g	Size: 58 x 38 x 20 cm Weight: kg
3660740013202	- pcs	3660740213206 -50 pcs
3660740013103	- pcs	3660740213107 -50 pcs
3660740013202	- pcs	3660740213206 -50 pcs

CUSTOM CODE: 90049010

## CLEANING &amp; STORAGE

**CLEANING:** The lenses should be cleaned regularly. Bollé Safety recommends only using the B-clean lens cleaning system.

**MAINTENANCE:** Pitted, scratched or damaged lenses must be replaced by replacing the complete product. If the frame is damaged also replace the complete product. After each use it is recommended to keep the spectacles or goggles in a microfiber pouch or rigid case.

**STORAGE:** All products must be stored in a dry and clean environment.



MADE IN: UK

DATASHEET: 01-07-2022

**bolle-safety.com**

**Europe Headquarters**  
34 rue de la Soie  
69100 Villeurbanne  
France  
Tél. : 00 33 (0) 4 78 85 23 64

**U.K. & Ireland**  
Unit C83 - Barwell Business Park  
Leatherhead Road - Chessington  
Surrey - KT9 2NY - UK  
Tél. : 00 44 (0)208 391 4700

**USA Office**  
2810 Caribou Court, Suite 160  
Carlsbad, CA 92010  
United States of America  
Tél. : 800 222 6553

**Asia-Pacific Office**  
32-40 Fairchild Street  
heatherton VIC 3202  
Australia  
Tél. : 1300 138 098



THE PPE EYEWEAR  
SPECIALISTS